Käyttöohje

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[Viitteet]

- Xu Xiaosi, Lv Shiming. The existence form of human chorionic gonadotropin in the body and its clinical application[J]. Foreign Medical Sciences (section of Clinical Biochemistry and Laboratory Medicine, 2005(11): 815-818.
- Wang Xiuping, Jiao Lin, Guo Hongyan. Progress in clinical application of human chorionic gonadotropin[J]. Pharmaceutical Care and Research, 2010, 10(03): 185-189.
- Zhang Ning, Lu Haiou, Ding Xiaoping. Potential sources of HCG during pregnancy and outside of pregnancy[J]. Foreign Medical Sciences-Section of Maternal and Child Health. 2003, 2(14): 5.
- Wang Hong, Li Hanbin, Shan Wenna, Ma Wencheng. Analysis of early pregnancy test strips for diagnosis of early pregnancy [J]. Heilongjiang Medical Journal, 1996(06): 41.
- Chen Hui. Application of human chorionic gonadotropin detection in obstetrics and gynecology [J]. China Practical Medicine, 2014, 9(15): 112-113.
- Xu Li. Analysis of human chorionic gonadotropin levels in early pregnancy[J]. China Continuing Medical Education, 2015, 7(14): 36-37.
- Wang Mei, Kang Hong, Li Yue. Exploration on the quality control method of colloidal gold immunochromatographic detection of urine hCG[J]. Medical Journal of West China, 2005(06): 105-106.
- Lee T. Human Chorionic Gonadotropin: Properties and Assay Methods. Semi Onool, 1995, 22:121-129.
- Blithe DL, Richards RG, Skarulis MC, Free alpha molecules from pregnancy stimulate secretion of prolactin from human decidual cells: a novel function for free in pregnancy. Endocrinology, 1991, 129:2257-2259.

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Jakelija Oriola Finland Oy Orionintie 5 02200 Espoo, Suomi www.pregcheck.fi

Guangzhou Decheng Biotechnology Co., Ltd. Floor 3/4/5/7, Building A1, No. 12, Nanyun 1st Road, Science City, Huangpu District, Guangzhou, Guangdong, 510663, P.R. Kiina www.dochebio.com service@dochebio.com

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[Referenser]

- [1] Xu Xiaosi, Lv Shiming. The existence form of human chorionic gonadotropin in the body and its clinical application[J]. Foreign Medical Sciences [section of Clinical Biochemistry and Laboratory Medicine, 2005(11): 815-818.
- [2] Wang Xiuping, Jiao Lin, Guo Hongyan. Progress in clinical application of human chorionic gonadotropin[J]. Pharmaceutical Care and Research, 2010, 10(03): 185-189.
- [3] Zhang Ning, Lu Haiou, Ding Xiaoping. Potential sources of HCG during pregnancy and outside of pregnancy[J]. Foreign Medical Sciences-Section of Maternal and Child Health. 2003, 2(14): 5.
- [4] Wang Hong, Li Hanbin, Shan Wenna, Ma Wencheng. Analysis of early pregnancy test strips for diagnosis of early pregnancy [J]. Heilongjiang Medical Journal, 1996(06): 41.
- [5] Chen Hui. Application of human chorionic gonadotropin detection in obstetrics and gynecology [J]. China Practical Medicine, 2014, 9(15): 112-113.
- [6] Xu Li. Analysis of human chorionic gonadotropin levels in early pregnancy[J]. China Continuing Medical Education, 2015, 7(14): 36-37.
- [7] Wang Mei, Kang Hong, Li Yue. Exploration on the quality control method of colloidal gold immunochromatographic detection of urine hCG[J]. Medical Journal of West China, 2005(06): 105-106.
- [8] Lee T. Human Chorionic Gonadotropin: Properties and Assay Methods. Semi Onool, 1995, 22:121-129.
- [9] Blithe DL, Richards RG, Skarulis MC, Free alpha molecules from pregnancy stimulate secretion of prolactin from human decidua cells: a novel function for free in pregnancy. Endocrinology, 1991, 129:2257-2259.

Importör ValMed A/S Bredgade 30 1260 København K, Denmark info@valmed.dk	Distributör Oriola Finland Oy Orionintie 5 02200 Espoo, Finland www.pregcheck.fi
--	--

	Guangzhou Decheng Biotechnology Co., Ltd. Floor 3/4/5/7, Building A1, No. 12, Nanyun 1st Road, Science City, Huangpu District, Guangzhou, Guangdong, 510663, P.R. China www.dochebio.com service@dochebio.com
---	--

	Qarad EC-REP BV Pas 257, 2440 Geel, Belgium
---	--

Instruction for Digital Pregnancy Rapid Test

INTENDED USE

This test is used for *in vitro* qualitative detection of Human Chorionic Gonadotropin (HCG) in human urine. This device is a self-testing IVD medical device for home use. It is suitable for the rapid diagnosis of pregnancy.

TESTING PRINCIPLE

Human Chorionic Gonadotropin (HCG) is a glycoprotein produced by the placenta during pregnancy. It can be detected in the urine of pregnant women. This test uses lateral flow immunoassay and light reflection for the detection of HCG in urine specimens. When the specimen is applied to the test, the HCG will react with the anti-β-HCG antibody-colloid gold conjugate and form a compound. As the liquid flows to the test area, the compound will be captured by the anti-α-HCG antibody immobilized on the test area, and a coloured line will form in the test area. The test will detect the light intensity by using an LED as the light source. After that, the result will be displayed on the display screen. The test diagnoses pregnancy quickly and easily. The test has a sensitivity of 25 mIU/ml. It can detect pregnancy as early as 1 day before the expected period.

COMPONENTS OF THE TEST KIT

REF 0811S7X001
1 x Test stick
1 x Instructions for Use
1 x Desiccant

Optional materials (not provided): A urine collection cup to collect the urine sample.

STORAGE CONDITIONS AND SHELF LIFE

1. Store in a dry place at 2°C - 30°C, do not freeze.
2. The shelf life is 30 months. You should use the test within one hour after opening the foil pouch.
3. See the product package for the manufacture date and expiration date.

SAMPLE REQUIREMENTS

You can either collect a urine sample into a urine collection cup or urinate directly on the absorbent tip of the test stick. Urine from any time of the day can be used, but morning urine is best. You should carry out the test immediately after urine collection. If you cannot carry out the test immediately, the urine can be stored at 2°C - 8°C for a maximum of 2 days. Make sure that the urine is at room temperature when applied to the test.

TEST METHOD

The test should be carried out at room temperature.


1. Take out the test stick from the foil pouch and use it within 1 hour of opening the foil pouch.
2. Testing from the urine stream: place the absorbent tip pointing downwards into your urine stream (Fig. 2). Or testing from collected urine: after collecting the urine sample, dip more than half of the absorbent tip into the urine (Fig. 3). Caution: In both cases, make sure that the rest of the test stick does not get wet.
3. After the urine has been applied on the absorbent tip, the hourglass symbol  will start flashing on the display screen (Fig. 4) to show that the test is working.
4. When you hear two beeps, remove the test stick from the urine and replace the cap, and place the test stick in a horizontal position. Caution: Never hold the test stick with the absorbent tip pointing upwards (Fig. 5).
5. Within 5 minutes, two beep sounds will be heard and your result will be shown on the display screen.

Fig.1 Components of the test stick


Fig.2 Test in urine stream

Remove the test stick from the urine stream when you hear two beeps.


Fig.4 Hourglass symbol

Or Fig.3 Test in a urine cup



Fig.5 Never hold the test stick with the absorbent tip pointing upwards

INTERPRETATION OF TEST RESULTS

Positive: If the screen displays "Pregnant" (positive), it indicates that you may be pregnant. Negative: If the screen displays "Not Pregnant" (negative), it indicates that you are not pregnant. Error: If the screen displays "Error" (invalid), it indicates that the test may be incorrect or invalid. Please check whether the test was performed in accordance with the instructions for use and repeat with a new test.

Pregnant

Positive

Not Pregnant

Negative

Error

Error

LIMITATIONS OF TEST METHOD

1. The product is only intended for HCG testing in human urine.
2. Unfresh or unclean test sample, or test sample with impurities, will produce incorrect results. Unclean urine collection containers, or improper operation will also produce incorrect results.
3. The test is intended for qualitative detection of pregnancy. A quantitative result should be further confirmed by other methods.
4. In some cases, the test results will be affected. For example: Injection or taking of hormone or steroid medications, those with history of previous injection of HCG, germ cell tumor of ovary (GCTO) and others.
5. The test results of this product will not be affected by haematuria, bilirubin or turbid urine. An elevated urine chyle level may affect the test result.
6. A false negative result may occur if the HCG level is too high, due to limitations of the test methodology.

PERFORMANCE CHARACTERISTICS

Sensitivity: 25 mIU/ml.

Specificity: This product does not produce a cross-reaction with luteinizing hormone (LH) (≤500 mIU/ml), follicle-stimulating hormone (FSH) (≤1000 mIU/ml) or thyroid-stimulating hormone (TSH) (≤1000 mIU/L).

PRECAUTIONS

1. This product is only intended for *in vitro* testing of human urine.
2. The product is for single use. It should be used before its expiry date. Reuse is not possible.
3. The test should be used indoors.
4. Please read the instructions for use carefully before use and follow the instructions, otherwise the test results will be affected.
5. The desiccant in the foil pouch is for storage purposes only. Do not swallow.
6. Do not use the test if the pouch is punctured or not well sealed.
7. When testing, make sure that urine is only applied to the absorbent tip to avoid causing false negative results.
8. If an "Error" symbol appears on the display screen, it indicates that the test is incorrect or invalid. Please make sure that the correct test method is followed and retest with a new test stick.
9. An hourglass symbol will appear on the display screen after the urine sample is applied. If it does not appear on the screen within one minute of adding the urine sample, it indicates that the operation of the test is abnormal. It is necessary to retest with a new test stick.
10. Urine collected at any time of the day can be used for the test. However, the first morning urine is the best as it contains the highest level of HCG and will increase the accuracy of the pregnancy test result.
11. Common medications (such as cold medicines, antibiotics, painkillers, etc.) do not affect the accuracy of the test results, but injecting or taking fertility drugs will affect the results.
12. A positive test result may indicate pregnancy. Please contact your healthcare provider to confirm the pregnancy and receive further advice.
13. Generally, a positive result indicates that you are pregnant. However, a positive result may also occur in the following situations: 1. An abnormal pregnancy, such as ectopic pregnancy, inevitable abortion, missed abortion or incomplete abortion. 2. Uterine tumours, hydatidiform moles or elevated levels of HCG in the urine related to menopause. To confirm pregnancy, further testing should be done by a doctor based on the user's medical history.
14. When the test result is negative, it indicates that you may not be pregnant. If your expected period has not come within one week, please test again. Please contact your healthcare provider if your second test is still negative.
15. If the test result is negative but pregnancy is suspected, retest with morning urine after 48 to 72 hours.
16. A hook effect (at 2-3 months of pregnancy) may occur when the HCG level is too high, causing a false-negative test result.
17. Before disposing of the product, remove the battery from the test stick. Insert a coin into the end of the stick and twist the cover open (Fig. 1). Dispose of the battery according to the local recycling protocol. Recycle the instructions for use as paper. Recycle the carton as cardboard. Warnings: Do not dispose of the battery via household waste. Do not disassemble or charge the battery. Do not destroy the battery by burning.
18. For any questions or further information about the use of this product, please contact the distributor.

FAULT & FAILURE ANALYSIS

Failure	Fault Analysis	Treatment
Display shows nothing.	No battery charge or the battery is in bad contact.	Please check the expiry date.
Error symbol indicating that the test is invalid.	The sample volume is too low or the instructions for use were not followed correctly.	Test again with a new test stick, following the instructions carefully.
The buzzer does not beep.	Bad contact in connection wire.	Contact the distributor.

KEY TO SYMBOLS USED

 Consult instructions for use	 Store at 2°C - 30°C
 Use-by date	 Manufacturer
 Batch code	 <i>In vitro</i> diagnostic medical device
 Do not use if package is damaged	 Authorized representative in the European Community
 Keep dry	 Keep away from sunlight
 Contains sufficient for <n> tests	 Catalogue number
 Do not reuse	 WEEE (Waste Electrical and Electronic Equipment)

[References]

- [1] Xu Xiaosi, Lv Shiming. The existence form of human chorionic gonadotropin in the body and its clinical application[J]. Foreign Medical Sciences [section of Clinical Biochemistry and Laboratory Medicine, 2005(11): 815-818.
- [2] Wang Xiuping, Jiao Lin, Guo Hongyan. Progress in clinical application of human chorionic gonadotropin[J]. Pharmaceutical Care and Research, 2010, 10(03): 185-189.
- [3] Zhang Ning, Lu Haiou, Ding Xiaoping. Potential sources of HCG during pregnancy and outside of pregnancy[J]. Foreign Medical Sciences-Section of Maternal and Child Health. 2003, 2(14): 5.
- [4] Wang Hong, Li Hanbin, Shan Wenna, Ma Wencheng. Analysis of early pregnancy test strips for diagnosis of early pregnancy [J]. Heilongjiang Medical Journal, 1996(06): 41.
- [5] Chen Hui. Application of human chorionic gonadotropin detection in obstetrics and gynecology [J]. China Practical Medicine, 2014, 9(15): 112-113.
- [6] Xu Li. Analysis of human chorionic gonadotropin levels in early pregnancy[J]. China Continuing Medical Education, 2015, 7(14): 36-37.
- [7] Wang Mei, Kang Hong, Li Yue. Exploration on the quality control method of colloidal gold immunochromatographic detection of urine hCG[J]. Medical Journal of West China, 2005(06): 105-106.
- [8] Lee T. Human Chorionic Gonadotropin: Properties and Assay Methods. Semi Onool, 1995, 22:121-129.
- [9] Blithe DL, Richards RG, Skarulis MC, Free alpha molecules from pregnancy stimulate secretion of prolactin from human decidua cells: a novel function for free in pregnancy. Endocrinology, 1991, 129:2257-2259.

Importer ValMed A/S Bredgade 30 1260 København K, Denmark info@valmed.dk	Distributör Oriola Finland Oy Orionintie 5 02200 Espoo, Finland www.pregcheck.fi
--	--

	Guangzhou Decheng Biotechnology Co., Ltd. Floor 3/4/5/7, Building A1, No. 12, Nanyun 1st Road, Science City, Huangpu District, Guangzhou, Guangdong, 510663, P.R. China www.dochebio.com service@dochebio.com
---	--

	Qarad EC-REP BV Pas 257, 2440 Geel, Belgium
---	--